



Care Cover Fabric Protection Shield

500mL Spray Bottle

The Care Cover Fabric Protection Shield, once applied, will assist in repelling stains on fabric furniture and surfaces, giving you time to remove marks and soiling of foreign substances before they penetrate the fibres of your fabric.

The Care Cover Fabric Protection Shield will be effective for three years from date of application (based on normal use of furniture such as couch or sofa).

The Care Cover Fabric Protection Shield will last through 3 professional cleans and up to three years, but can be re-applied as required by following the steps below.

Directions for use:

- 1 - Clean the fabric, ensuring no marks, dust or residue remain before applying Care Cover Protection Shield.
- 2 - Moving left to right and from at least 50cm away, spray an even mist over the surface of the fabric.
- 3 - Allow mist to absorb into the fibres, ensuring high use areas of the fabric are well covered.
- 4 - Repeat steps 2 and 3 until all fabric has been completely covered.
- 5 - Allow to cure for 24 hours (longer in colder months/climates).

Includes:

- 500mL spray bottle

Application of the Protection Shield has proven to help prevent liquid from absorbing into and staining fabric.

Opposite image shows reaction of liquids such as coffee and water when poured onto a treated fabric surface.

